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E-mail: enquiries@handicare.co.uk

Website: www.handicare.co.uk

Important information about your Minivator Bath Lift

Figure 1

As the owner of a Minivator bath lift, (Figure 1) we would like to make you aware of a rare potential fault which has occurred on a small number of these devices whilst in use.

There have been a very small number of reported faults with units using the style of handset ('keyhole' shape, blue or grey in colour) shown in **Figure 2**, where the bath lift would not rise-up out of the bath when the handset up button was pressed.

Although this problem is very rare, Handicare considers your safety to be paramount and we are therefore taking the following actions to ensure that you can continue to bathe with confidence when using your bath lift.

You will be contacted in the near future so we can help you exchange your existing bath lift handset for a new one (**Figure 3**). You will get a new handset **totally free of charge**, we will make all the necessary arrangements to have the handset replaced; you do not need to do anything at this time.

Prior to receiving your new handset we would ask you to bathe safely by following the advice below:

- When using the bath lift, please ensure that you have somebody else in the house that can be called to assist should you have any problems.
- If you live alone, please ensure that you have the means to call for assistance whilst in the bath should this be required. This could be a warden call system, telephone, mobile phone or other communication device.
- The fault has been identified to water leaking into a small number of handsets during normal use. To minimise the risk of this occurring prior to your handset being exchanged, please try and keep the bathlift handset out of the bath water during use. By attaching the handset to the top rim of the bath using the fitted suction cups you will achieve this.

As we have indicated, it is very rare for a fault to occur and we expect your bath lift to continue to provide you with safe and trouble free use. The replacement of your bath lift handset is a precaution to ensure we completely remove the potential of this occurring.

Figure 2



Figure 3

Note: If your Bath lift already has the handset shown in Figure 3 you do not need to take any action.

A5 first class postage paid Bath Lift Status Card

Our records indicate that a Minivator Bath Bliss Bath Lift was delivered to this address. This unit may require a safety upgrade. We would be grateful if you could complete and return this postage paid card as appropriate to confirm whether or not you have the unit.

If the Bath Lift is still at this address please write the serial number in the box provided below.

The serial number is found on label similar to this.

The label is on the back of bath lift not the handset.

If the Bath Lift is in use at another address please provide contact details here if you have them:

Name of user _____

Address _____

Postcode _____

Tel _____

If the Bath lift is no longer at this address and you are unaware of its whereabouts please tick this box ☐

Once completed please return this card as soon as possible.

Thank you for your assistance. If you have any questions regarding this matter or a Minivator Bath Lift please contact Handicare 0844 225 3121.

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XX July, 2012

Dear Valued Partner,

This letter is to inform you of an amendment to the field safety corrective action on the Legendair™ and Supportair™ ventilators, initiated by Covidien in March 2012, related to a potential failure of a component identified as capacitor C53.

This occurrence may result in the inability of the ventilator to switch automatically from A/C mains power to internal battery operation when the ventilator is disconnected from A/C mains during use or when there is a failure of the AC mains power during use. The ventilator will run on external battery or A/C mains power.

The affected ventilator codes are referenced below:

| Legendair™ Single Branch | | Legendair™ Double Branch | | Supportair™ | |
|--------------------------|-------------------------------|--------------------------|-------------------------------|-------------|-----------------|
| 4095700 | INTL LEGENDAIR SINGLE BRANCH | 4095900 | INTL LEGENDAIR DOUBLE BRANCH | 4096200 | INTL SUPPORTAIR |
| 4095701 | GB NO LEGENDAIR SINGLE BRANCH | 4095901 | GB NO LEGENDAIR DOUBLE BRANCH | 4096201 | GB SUPPORTAIR |
| 4095702 | DE LEGENDAIR SINGLE BRANCH | 4095902 | DE LEGENDAIR DOUBLE BRANCH | 4096202 | DE SUPPORTAIR |
| 4095705 | CA LEGENDAIR SINGLE BRANCH | 4095905 | CA LEGENDAIR DOUBLE BRANCH | | |
| 4095706 | JP LEGENDAIR SINGLE BRANCH | 4095906 | JP LEGENDAIR DOUBLE BRANCH | | |
| 4095709 | DK LEGENDAIR SINGLE BRANCH | 4095909 | DK LEGENDAIR DOUBLE BRANCH | | |
| 4095721 | AU LEGENDAIR SINGLE BRANCH | 4095921 | AU LEGENDAIR DOUBLE BRANCH | | |

According to our records, you have a ventilator that may be impacted by this issue.

Working in conjunction with several regulatory agencies, we have decided to further segment the population involved in this action into ventilator dependent and non-ventilator dependent patients. For ventilator dependent patients*, Covidien recommends replacement of the PCB (motor control board) which contains the C53 capacitor as soon as possible. Until the PCB is replaced, those patients who are equipped with the OpenPack™ or PowerPack™ external battery should keep the external battery connected at all times. Additionally, if a second ventilator has been provided per prescription, it should also be ready and available.

Attached to this field safety notice is the letter (Attachment A) we request you to send to all your end-users to inform them of the issue and to reassure them that by working with you, they can maintain proper operation of their ventilator until the PCB is replaced.

Because PCBs need to stay available for ventilator dependent patients, the replacement of the C53 capacitor for non-ventilator dependent patients needs to be included as part of the *revised* annual preventive maintenance. Replacement of the C53 component should occur no later than the next scheduled maintenance and should then be replaced annually as part of the

preventative maintenance schedule if PCBs are not available. Replacement of the PCB is the preferred action.

Instructions for replacing both the motor control board and the capacitor are provided in Attachment B. We request you to complete attached Verification Form provided in Attachment C and complete this with the serial numbers you have refurbished. If you have purchased a large number of devices, please complete this form regularly as you perform maintenance for your device and send the form back to the contact mentioned on the Verification Form.

If you do not perform your own preventive maintenance, please contact your maintenance provider to inform them of this update, or contact your local Covidien Service Center.

This notice is being distributed with the knowledge of appropriate regulatory authorities.


This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please do not hesitate to contact your Covidien representative or Technical Service if you have any additional questions.

Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.

We sincerely apologize for any inconvenience this may cause and appreciate your prompt attention to this matter.

Sincerely,

Michael A. Ronnigen
VP, Regulatory Affairs
Respiratory and Monitoring Solutions
Covidien

*Ventilator dependent patient = patient ventilated 16 or more hours/day

Non-ventilator dependent patient = patient ventilated less than 16 hours/day

Attachment A: End-User letter

Attachment B: Technical Service Bulletin

Attachment C: Verification form

Dear Valued Customer,

This letter is to inform you that Covidien, the manufacturer of the Legendair™ ventilator, has discovered that one of the components located on the Legendair™ and Supportair™ electronic Motor Control Board might exhibit degradation over time. The specific affected component is the C53 Capacitor. No failures have occurred due to component degradation but Covidien has decided to initiate action to prevent any failures from degradation over time. There are two actions that may be taken to ensure ventilators continue to function. Only one of these actions is needed address the issue:

1. Replace the affected motor board by a new board; or
2. Replace the affected C53 component yearly as part of preventive maintenance.

Your homecare provider will contact you to organize this modification.

The function of the C53 Capacitor is to switch the ventilator operation from main electric power to battery operation when the electrical plug is removed from the electrical wall outlet.

If the capacitor fails due to degradation, the switch from main electric power to battery operation might not occur, resulting in the ventilator ceasing to function. In this event an audible, non-resettable audible alarm will sound for a minimum of 10 seconds.

For this reason Covidien recommends following the precautions until your home care provider exchanges the board or replaces the component as part of preventive maintenance.

- 1) Make sure that the ventilator's internal battery is fully charged and functional before the ventilator is used, even if the ventilator will be plugged into the wall. If you unplug the ventilator from the wall outlet, please assess and confirm the continuation of the ventilator operation.
- 2) If you own an OpenPack™ or PowerPack™ external battery pack, leave the ventilator connected to the external battery pack at all times. Those battery packs will manage the ventilator power supply. See User's Manual for operating parameters.
- 3) Even when the patient is not connected to the ventilator, it is recommended to keep the ventilator plugged to the battery pack and into the wall outlet, so that the external battery is being continually recharged. Also, if the ventilator has been stored for a period of time, it is recommended that you recharge the internal battery before the ventilator is used.
- 4) If you own a second ventilator, make sure it is available and ready for use. It is recommended to keep the ventilator plugged into the wall outlet, so that the internal battery is being recharged. Also, if the ventilator has been stored for a period of time, it is recommended that you recharge the internal battery before the ventilator is used.
- 5) The Legendair ventilator continuously and automatically controls the state of the internal battery, even if the ventilator is being operated on wall power. However, it is recommended that you check the charge status of the internal battery once each month. This is a simple procedure that is described in the "Internal Battery Maintenance" section of the User's Manual.

Should you have any questions regarding this letter, please contact our office at [phone number]. .

Sincerely,

Home Healthcare Provider

**TSB1000014-01 Rev. B
ECO-R206600****International Technical Service Bulletin, Portable Ventilators,
Legendair and Supportair
Revised -Addendum to the Fields Safety Corrective Action**

This Technical Service Bulletin details the different actions that may be necessary to take in relation to the Urgent Field Safety Notice dd June 2012.

The affected devices cover the following legacy ranges:

Airox Legendair™ Single Branch
Airox Legendair™ Double Branch
Airox Supportair™


The change is focused on the Turbine Controller Printed Circuit Board fitted within these devices and specifically on the replacement of capacitor C53.

The following options will aid you in replacing the C53 or updating your device:


1. The New PCB 3832799 can be installed as per the service manual
 - Legendair Service Manual Rev B. Chapter 5 page 11
 - Supportair Service Manual Rev A. Chapter 5 page 12
2. Replacement of the existing 3812799 PCB that contains a new C53 component, is done following the above mentioned service manual references.
3. Replacement of the C53 capacitor on the existing PCB can be performed by the HCP or at a Covidien Technical Services location.

The capacitor has been released as a spare part to allow replacement of this capacitor, part number 10049759, which will be available for order in a pack containing 10 capacitors.

This change is limited to the devices with either Turbine Controller PCB 3812700 or 3812799 installed, and once the device has the new version Turbine Controller PCB 3832799 installed, no further action is required.


This PCB needs C53 replaced

This PCB requires no action

With the top cover removed the Turbine Controller PCB is located here.

The C53 is located here

This matrix shows which device and printed circuit board combination is impacted by this improvement.

If your device shows "Yes" you will need to replace the C53 as detailed below.

If your device has "No" within the matrix no action is required.

| Device | Turbine Controller PCB Installed | | |
|------------|----------------------------------|---------|---------|
| | 3812700 | 3812799 | 3832799 |
| Legendair | Yes | Yes | No |
| Supportair | Yes | Yes | No |

The updated preventative maintenance schedule is detailed below, only applies to devices with 3812700 or 3812799 installed.

This schedule replaces the one detailed in Chapter 2 of both the Legendair revision B and Supportair revision A Service Manuals.

Legendair

| Interval | Prefix |
|---------------|--------|
| Every Year | V3 |
| Every 3 years | V4 |
| Every 6 years | V5 |

| Operation | V3 | V4 | V5 |
|--|----|----|----|
| Checking user/machine interfaces, adjustment and calibration | X | | |
| Intentional vent stop alarm test | X | | |
| Flow laminator cleaning | X | | |
| Replacement of C53 capacitor | X | | |
| Replacement of internal battery | | X | |
| Replacement of turbine box | | | X |

Supportair

| Interval | Prefix |
|---------------|--------|
| Every Year | V2 |
| Every 3 years | V3 |
| Every 6 years | V4 |

| Operation | V2 | V3 | V4 |
|--|----|----|----|
| Checking user/machine interfaces, adjustment and calibration | X | | |
| Inspection date of oxygen hose | X | | |
| Replacement of C53 capacitor | X | | |
| Replacement of internal battery | | X | |
| Flow Laminator cleaning | | X | |
| Checking, adjustment and calibration O2 | | X | |
| Replacement of pilot solenoid valve | | X | |
| Replacement of complete O2 block | | | X |
| Replacement of cooling fan | | | X |
| Replacement of turbine box | | | X |

This change should be made with immediate effect.

Adam Pearce
EMEA Technical Product Manager
Covidien

URGENT MEDICAL DEVICE FSCA VERIFICATION FORM

| Customer Contact Details | Covidien Contact Details |
|---|--|
| Hospital / HCP: Covidien Account Number: | To: [please insert name Covidien commercial office] |
| Address: | Address: [please insert Covidien address] |
| Telephone n°: | Telephone n°: [please insert Covidien telephone number] |
| Fax n°: | Fax n°: [please insert Covidien fax number] |
| E-mail: | E-mail: [please insert contact e-mail address] |

I acknowledge receipt of the Field Safety Notice regarding Legendair™ and Supportair™ C53 FSCA _____ (signature)

Please complete below column, indicating all serial numbers you have serviced and indicate whether the C53 capacitor or Motor Control Board (MCB) was replaced.

| Serial numbers Legendair™ or Supportair™ ventilator | C53 component replaced on board (P/N 100049759) Yes/No | MCB with new C53 (P/N 3812799) Yes/no | New revision MCB (P/N 3832799) Yes/No |
|---|---|---|---|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

Please fax this form to the fax number referenced at the top of this form.

Name:
(please print)

Signature:

Date:

| FSN n° 2014-13659 | |
|-------------------|---------------------------|
| File | 13659 |
| Manufacturer | CircuLite GmbH |
| Device | Synergy Micro-pump system |
| Commission date | 08/05/2014 |

| | |
|------------------|---|
| FSN n°2014_13854 | |
| File | 13854 |
| Manufacturer | Becton Dickinson |
| Device | BD Epilor 7ml Luer Lok Loss of Resistance Syringe |
| Commission date | 08-05-2014 |

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FSN only available in English.

**To the ATTENTION of:
Operating Room Manager**

25 September 2013

URGENT MEDICAL DEVICE PRODUCT REMOVAL

Part Description / Part Number

| Part Number | Part Description | Lot Number |
|----------------|---|-------------------------------|
| 05.001.178.01S | Irrigation Tube Set, sterile, single pack | 3755327 7823268 7823269 |
| 05.001.178.05S | Irrigation Tube Set, sterile, pack of 5 units | 3755329 7884993 |

Dear Madam / Dear Sir,

Synthes is initiating a voluntary recall of the above mentioned articles and lots of the Irrigation Tube Set, sterile, due to the possibility of a leakage. Our records indicate that you may have inventory that is impacted by this removal.

Description of problem:

Synthes received a customer complaint for one device reporting a leak. No patient harm was reported.

This leak occurred at the coupling that connects the smaller and larger irrigation tubes, approximately 35cm from the working end of the tubing. The leak was due to an insufficient amount of adhesive applied in the manufacturing process to secure the smaller tubing into the connector. Synthes performed internal testing and confirmed that this issue affected other units in the above mentioned lots.

User risk:

Defective tubing not detected prior to use may be introduced for use with the drill into the operative field. In a worst case scenario, there is the potential for the tubing to leak into the operative field and mix with body fluids. With enough pressure, there is the potential for those fluids to splash back onto the user. If the user is not wearing appropriate PPEs (Personal Protective Equipment), like eye shields, then fluid exposure may occur. If exposure occurs the user would then require immediate irrigation of the exposed area and potential additional

testing and monitoring. Additional medical treatment based upon patient specific findings and the institution's biohazard exposure plan may be required.

Customer immediate actions:

1. Immediately identify and quarantine all unused products listed above in a manner that ensures the affected products will not be used.
2. Review, complete, sign and return the attached reply form to your local Synthes sales organisation in accordance with the directions on the form within 5 business days of receipt of this notification.
3. Return any affected product within 30 business days. A credit note will be issued for the returned items.
4. Forward this notice to anyone in your facility that needs to be informed.
5. If any product listed below has been forwarded to another facility, contact that facility to arrange return.
6. Maintain awareness of this notice until all products listed below have been returned to Synthes GmbH.
7. Maintain a copy of this notice.

The applicable regulatory agencies are being notified. Synthes GmbH is voluntarily taking this action.

We apologise for any inconvenience that this product removal may create and appreciate your cooperation with our request. Should you have any queries please do not hesitate to contact your Synthes sales representative.

Page 3 of 4

Thank you for your attention and cooperation.

Synthes GmbH

Claudia Allemann
Field Action Manager

Markus Wien
Director Quality Assurance Operations

Cc:

NOTICE: MEDICAL DEVICE REMOVAL R2013563**Verification Section**

| Part Number | Part Description | Lot Number |
|----------------|---|-------------------------------|
| 05.001.178.01S | Irrigation Tube Set, sterile, single pack | 3755327 7823268 7823269 |
| 05.001.178.05S | Irrigation Tube Set, sterile, pack of 5 units | 3755329 7884993 |

- ☐ We have located the identified product in stock; returned quantity is documented below, and a copy of this letter is being retained for our records.
- ☐ We do not have any identified product in stock; returned quantity is zero. We have retained a copy of this letter for our records.

RETURNED DEVICES (including quantity) and/or COMMENTS:

Hospital name: _____

Name/Title (please print) _____


Phone Number: _____

Signature and Date: _____

| FSN n° 2014-14482 | |
|-------------------|------------|
| File | 14482 |
| Manufacturer | Elekta Ltd |
| Device | RT Desktop |
| Commission date | 08/05/2014 |

Important Field Safety Notice

Recommended quality assurance after a calibration procedure

Product: Desktop Pro™

Date: November 2013

FCO: 200 01 501 041

ELEKTA

Important Field Safety Notice

We are providing the information in this Notice to notify you of an important safety issue that may exist on your equipment, and to inform you of any actions needed to safeguard both your staff and your patients. We ask that you please read and understand the content of this notice and implement any recommendations provided.

We also need you to acknowledge and accept this FCO by signing and returning the statement on the FCO Action Notification Report.

We advise you to insert this Notice in the applicable copy of the User Manual.

Recommended quality assurance after a calibration procedure

Product: Desktop Pro™

Date: November 2013

FCO: 200 01 501 041

| | |
|------------------|---|
| Scope: | Desktop Pro™ R4.2, R5.0, R5.1, R6.0, R6.1 |
| Description: | <p>During the treatment table, gantry, and collimator calibration procedures, it is possible to type incorrect values. If you do not do the quality assurance (QA) procedures after the calibration, it is possible that you will not find the error made during the calibration. For example:</p> <ol style="list-style-type: none"> 1. The Service User types the minimum and maximum lateral values for the treatment table in the incorrect sequence. 2. You do not do the post-calibration QA checks, or you do the QA checks incorrectly. 3. You use the machine for clinical operation. 4. You do not do the daily QA checks, or you do the QA checks incorrectly. 5. When you do a remote automatic table movement during an online XVI workflow, you do not see that the treatment table moves in the opposite lateral direction. 6. You deliver the treatment with the patient in the incorrect position. |
| Clinical impact: | If the treatment table, gantry, or collimator calibration is incorrect, and you do not do the recommended QA procedures, it can cause clinical mistreatment. |

This Notice has been notified to the appropriate Regulatory Authorities

Important Field Safety Notice

| | |
|----------------------|---|
| Solution: | <p>Elekta recommends that:</p> <ul style="list-style-type: none"> You do the daily QA procedures in the <i>Precise Table Instructions for Use (IFU) and Desktop Pro™ IFU</i>. You do a daily QA check to make sure that the magnitude and direction of the treatment table movements are correct. If you have XVI, you do the QA procedures in the <i>XVI IFU</i> and the procedure in this Notice to do the daily phantom scan. |
| Technical Reference: | CLM 01681502 |
| Contact: | If you have any queries about this Notice, please contact your local Elekta office. |

Recommended QA for XVI

The XVI system is a precision instrument, and it is necessary to do regular planned maintenance, user maintenance, safety checks, and QA checks. These maintenance tasks and checks help to give safe and satisfactory system operation with the necessary reliability.

It is the responsibility of the Clinical User to make sure that the QA checks are completed before clinical operation.

It is the responsibility of the authority that has the control of the equipment to:

- Set up a QA program to make sure that the image quality is satisfactory, the equipment aligns correctly, and the calibration is correct.
- Make sure that the QA program obeys applicable laws and regulations in the jurisdiction of the equipment.
- Make sure that the QA checks are completed.
- Keep the necessary QA records.

This Notice has been notified to the appropriate Regulatory Authorities

Important Field Safety Notice

Elekta recommends that the authority that has the control of the equipment uses these documents to help set up a QA program:

- Klein et al., 2009. 'Task Group 142 report: Quality assurance of medical accelerators.' Medical Physics, 36(nine), pp. 4197-4212.
- Bissonnette et al., 2008. 'Quality assurance for the geometric accuracy of cone-beam CT guidance in radiation therapy.' Int. J. Radiat. Oncol. Biol. Phys., 71(1), pp.S57-S61.
- Bissonnette et al., 2008. 'A quality assurance program for image quality of cone-beam CT guidance in radiation therapy.' Med. Phys., 35(5), pp.1807-1815.
- Sykes et al., 2008. 'Measurement of cone beam CT coincidence with megavoltage isocentre and image sharpness using the QUASAR™ Penta-Guide phantom.' Phys. Med. Biol., 53, pp.5275-5293.
- Cho et al., 2005. 'Accurate technique for complete geometric calibration of cone-beam computed tomography systems.' Med. Phys., 32(4), pp.968-983.

Doing the daily phantom scan

This is a procedure to make sure that the communication between the Table Move Assistant and the treatment table is active.

Elekta recommends that you do a daily phantom scan before the daily treatment cycle. This will make sure that communication between the Table Move Assistant and the treatment table is active for the radiation treatments.

1. Do the daily phantom scan, see the documents that follow for more information:
 - The AAPM TG-179 report, Bissonnette J-P et al, Quality assurance for image-guided radiation therapy utilizing CT-based technologies: A report of the AAPM TG-179. This is available at <http://www.aapm.org>.
 - Kurt Stump, Elekta XVI Volumetric Imaging and Registration QA With the MIMI Phantom. This is available at <http://www.standardimaging.com>.

This procedure also identifies an incorrect calibration of the treatment table.

This Notice has been notified to the appropriate Regulatory Authorities

Important Field Safety Notice

Safety reference

The following warnings and cautions are associated with this notice:

WARNING

Do not use the equipment unless the user maintenance, safety checks, and QA checks are satisfactorily completed, and the planned maintenance program is completed to the current date. If you ignore this warning, it can cause fatal injury, clinical mistreatment, and damage to the equipment.

This Notice has been notified to the appropriate Regulatory Authorities

FCO Action Notification Report

Please complete the details below and sign the appropriate acknowledgement section:

- Existing installations; Acknowledgement by the customer
- New installations: New installation confirmation by the installing Elekta or Representative employee

Please return this report to your local Elekta Office or Representative, as soon as possible and within 30 days at the latest.

***The information in this FCO has been provided to address a safety issue and therefore the customer is expected to acknowledge and accept the recommendations given, and ensure they are implemented. By refusing to implement the recommendations, the customer assumes full responsibility and liability for all matters (including costs, losses, claims, and expenses) resulting, whether directly or indirectly from not implementing such recommendations. Further the customer will hold Elekta harmless from all matters (including costs, losses, claims and expenses) resulting, whether directly or indirectly from not implementing such recommendations. Failure to sign and return the acknowledgement may affect any follow-up actions necessary for us to take, and may require Elekta to report to the Regulatory Authorities in your country.**

| | | | |
|------------------|---|----------|----------------|
| Classification: | Important Field Safety Notice | FCO Ref: | 200 01 501 041 |
| FCO description: | Recommended quality assurance after a calibration procedure | | |
| Scope: | Desktop Pro™ R4.2, R5.0, R5.1, R6.0, R6.1 | | |

| | |
|--|-------------------------|
| Hospital: | |
| Device Serial No: (e.g. linac - if applicable) | Location or Site No: |

| | |
|--|--------|
| Acknowledgement by customer*: This notification to be signed by the customer I acknowledge that I have read and understood this Notice and accept implementation of any given recommendations: | |
| Name: | Title: |
| Signature: | Date: |

| | |
|--|--------|
| New installation confirmation: This notification to be signed by the installing Elekta or Representative employee | |
| Name: | Title: |
| Signature: | Date: |

This Notice has been notified to the appropriate Regulatory Authorities

| FSN n° 2014-14501 | |
|-------------------|-------------|
| File | 14501 |
| Manufacturer | Siemens AG |
| Device | Artis zeego |
| Commission date | 08/05/2014 |

Siemens AG, H IM AX QM, Siemensstr. 1, 91301 Forchheim, Germany

To all users of Artis systems with software version Artis VC21B (without patch 2) configured with a Large Display

| | |
|------------|-----------------------------------|
| Name | Michael Greiffenhagen |
| Department | H IM AX MK IPM |
| Phone | +49 (9191) 18-6540 |
| Fax | +49 (9191) 18-9912 |
| Mobile | +49 (173) 2561023 |
| Email | michael.greiffenhagen@siemens.com |
| Date | 2013-10-25 |

Important customer safety notice

Update AX048/13/S

Safety-related information regarding corrective measures for Artis systems running software version VC21B in conjunction with a Large Display

Dear Customer,

This is to inform you about a possible problem with Artis systems that are running software version VC21B in conjunction with a Large Display.

When does this problem arise and what possible risks are there?

If the system is not turned "OFF" manually at regular intervals, an automatic reset of the Large Display Controller can occur. During the reset the system switches for approximately 1-2 minutes into the "Bypass Fluoro" mode. During this time, only the native fluoroscopic image is displayed.

What measures are to be taken to avoid possible risks?

The issue described above can be resolved by installing Update AX048/13/S, distributed with this notice.

Until Update AX048/13/S has been installed, we recommend that you proceed as follows (once a day for monoplane systems and after every patient examination for biplane systems): Turn the system "OFF" (**„System shutdown“**) and then turn the system "ON" again via the user interface (**power cycle**). Performing a "System restart" or "application restart" is not sufficient to prevent the problem from occurring.

Siemens AG
Healthcare Sector; Management: Hermann Requardt
Imaging & Therapy Division Systems; Management: Bernd Montag
Angiography & Interventional X-Ray Systems; Management: Heinrich Kolem

Siemensstr. 1
91301 Forchheim
Germany

Tel.: +49 (9191) 18 0
Fax: +49 (9191) 18 9999

Siemens Aktiengesellschaft: Chairman of the Supervisory Board: Gerhard Cromme;
Managing Board: Joe Kaeser, Chairman, President and Chief Executive Officer; Roland Busch, Klaus Helmrich,
Barbara Kux, Hermann Requardt, Siegfried Russwurm, Peter Y. Solmssen, Michael Suess, Ralf P. Thomas
Registered offices: Berlin and Munich, Germany; Commercial registries: Berlin Charlottenburg, HRB 12300, Munich, HRB 6684
WEEE-Reg.-No. DE 23691322

What risks are there for patients who have previously been examined or treated using this system?

We do not consider it necessary in this case to re-examine patients.

We thank you for your cooperation in dealing with this safety notice and request that you promptly notify and instruct accordingly all the staff at your organization who have to be aware of this problem. Please also forward this safety information to other organizations that could be affected by this measure. Please observe this safety notice, and comply with the corresponding measures until the update has been fully completed.


If the device has been sold and is therefore no longer in your possession, please forward this safety notice to the new owner. We would also request that you inform us of the identity of the device's new owner where possible.

Best regards,
SIEMENS AG Healthcare Sector
Business Unit AX

Dr. Heinrich Kolem
CEO HIMAX

Rudolf Glaumner

Wolfgang Hofmann
Safety Officer Medical Device

| FSN n° 2014-14522 | |
|------------------------|--|
| File | 14522 |
| Manufacturer | Edwards Lifesciences |
| Device | Peripheral Access Arterial / Femoral Arterial Cannulas (37695); Peripheral Access Venous / Femoral Venous Cannulas (34905); |
| Commission date | 08/05/2014 |

Urgent Field Safety Notice

ADEPT® Hip Resurfacing System

FSCA Identifier: FSN-004517

Type of Action: Field Safety Notice.

Date: 18 December 2013

Attention:

Trust Chief Executives, Clinical Director- Orthopaedics Department, Orthopaedics Theatre Manager, Safety Liaison Officer, General Managers – Private Sector Hospitals.

Details of affected devices:

ADEPT® Hip Resurfacing System for head sizes under 48mm.

| Catalogue Number | Description |
|-------------------------|--|
| 172-445 | ADEPT® Metal on Metal Cup 44mm with 38mm bore |
| 172-461 | ADEPT® Metal on Metal Cup 46mm with 38mm bore |
| 172-465 | ADEPT® Metal on Metal Cup 46mm with 40mm bore |
| 172-481 | ADEPT® Metal on Metal Cup 48mm with 40mm bore |
| 172-485 | ADEPT® Metal on Metal Cup 48mm with 42mm bore |
| 172-501 | ADEPT® Metal on Metal Cup 50mm with 42mm bore |
| 172-505 | ADEPT® Metal on Metal Cup 50mm with 44mm bore |
| 172-521 | ADEPT® Metal on Metal Cup 52mm with 44mm bore |
| 172-525 | ADEPT® Metal on Metal Cup 52mm with 46mm bore |
| 172-541 | ADEPT® Metal on Metal Cup 54mm with 46mm bore |
| 174-461 | ADEPT® Extra Fixation Metal on Metal Cup 46mm with 38mm bore |
| 174-501 | ADEPT® Extra Fixation Metal on Metal Cup 50mm with 42mm bore |
| 174-541 | ADEPT® Extra Fixation Metal on Metal Cup 54mm with 46mm bore |
| 181-138 | ADEPT® Metal on Metal Resurfacing Head 38mm |
| 181-140 | ADEPT® Metal on Metal Resurfacing Head 40mm |
| 181-142 | ADEPT® Metal on Metal Resurfacing Head 42mm |
| 181-144 | ADEPT® Metal on Metal Resurfacing Head 44mm |
| 181-146 | ADEPT® Metal on Metal Resurfacing Head 46mm |

Forever Active

MatOrtho Limited | 13 Mole Business Park | Randalls Road | Leatherhead | Surrey | KT22 7BA | United Kingdom

T: +44 (0)1372 224 200 | F: +44 (0)1372 366 343 | info@MatOrtho.com | www.MatOrtho.com

Company registered in England: 7323441

Description of Problem:

MatOrtho has added a warning to the IFU for the **ADEPT® Hip Resurfacing** System that the revision rates are greater for smaller head/cup sizes and has changed the indications to clarify that hip resurfacing is for younger or more active patients where conventional total hip replacement is likely to require revision.

A warning has also been added to the IFU to reflect that modular heads have not been available for this system since early 2012 and therefore if the resurfacing head must be revised to a total hip, the acetabular cup should also be revised, even if it is well fixed.

The NJR 10th Annual Report Table 3.11 shows a Kaplan-Meier estimate of the cumulative percentage probability of first revision (95% CI) for the **ADEPT® Hip Resurfacing** System of 8.7 (6.91-10.93) at 7 years. This is starting to trend slightly outside the NICE guidelines of 10% at 10 years.

Independent analysis commissioned by MatOrtho of the NJR data downloaded 29 August 2013 shows that for head sizes 44mm and above, the estimate of the revision rate at 7 years is less than 7.0%

In general, revision rates for smaller sizes are higher than for larger sizes although the data is weak due to low numbers.

We are aware that other manufacturers are recommending head sizes 48mm and above and without robust evidence to support acting differently, we have decided to only supply head sizes 48mm and above.

The data available for **ADEPT® Hip Resurfacing** System is not robust enough to provide positive or negative conclusions for various reasons including:

1. Indications and clinical practice for all metal on metal resurfacing devices has changed over the period that the data has been collected and does not reflect current use.
2. The overall performance of the metal on metal resurfacing as a class of various emerging good and indifferent devices is confused by the variable performance of the different devices available during the NJR study period. **ADEPT® Hip Resurfacing** System is one of the two best systems available according to various registries.
3. The dataset available for **ADEPT® Hip Resurfacing** System is small and therefore the confidence intervals are wide. The more specific the analysis, the less statistically robust it becomes. Larger sizes in younger patients are however tracking towards acceptably low revision rates.

MatOrtho has decided to clarify the **ADEPT® Hip Resurfacing** System indications for use to younger or more active patients and to warn about smaller implant sizes and to underline this, MatOrtho will voluntarily withdraw the above listed **ADEPT® Hip Resurfacing** System products.

Forever Active

MatOrtho will first inform users directly then exchange instrument sets so that cutters and trials for the withdrawn sizes will not be available and withdraw the above implant products. This will be done on a centre by centre basis and will be coordinated by the local MatOrtho representative worldwide.

Advice on Action to be taken by the user:

The above devices should not be used and should be quarantined, ready for collection by the MatOrtho representative.

Be aware that the revision rate for smaller head sizes is greater than for larger head sizes for both males and females.

All patients receiving metal on metal devices should be followed as per local orthopaedic association guidelines or the guidelines produced by MHRA.

Transmission of the Field Safety Notice:

This notice needs to be passed to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Contact reference person:

Andrew Woodall
QA and RA Director
MatOrtho Ltd
13 Mole Business Park
Randalls Road
Leatherhead
Surrey
KT22 7BA

Telephone: +44 (0)1372 366 300

Email: adept.update@matortho.com

Forever Active

The undersigned confirms that this notice has been notified to the appropriate Regulatory Agency.

The data currently available demonstrates that the **ADEPT® Hip Resurfacing** System continues to offer an effective treatment for appropriate patients when current guidelines are followed.

To confirm receipt of this FSN, please complete and return the acknowledgment in appendix A.

Signed:

Date:

QA and RA Director

Forever **Active**

MatOrtho Limited | 13 Mole Business Park | Randalls Road | Leatherhead | Surrey | KT22 7BA | United Kingdom

T: +44 (0)1372 224 200 | F: +44 (0)1372 366 343 | info@MatOrtho.com | www.MatOrtho.com

Company registered in England: 7323441

Appendix A

This letter acknowledges receipt of Field Safety Notice (FSN-004517) dated 18 December 2013 issue by MatOrtho Limited.

Please fax or email this completed document to MatOrtho at +44 (0)1372 366343 or adept.update@matortho.com

| | |
|-----------------------------------|-------|
| Signature | <hr/> |
| Print Name | <hr/> |
| Hospital | <hr/> |
| City | <hr/> |
| Country | <hr/> |
| Telephone Number or email address | <hr/> |

Forever Active

| FSN n° 2014-14639 | |
|-------------------|----------------------------|
| File | 14639 |
| Manufacturer | Palodex Group Oy |
| Device | FOCUS(TM) -Intraoral X-ray |
| Commission date | 08/05/2014 |

FIELD SAFETY NOTICE

Product: FOCUS™ 3.0 Intraoral X-ray unit (type: FOCUS)
Affected devices: F21221 – F21643 (excluding F21642)

FOCUS™ 3.0 X-ray unit field inspection required

We regret to inform you that Instrumentarium Dental has delivered a number of FOCUS™ 3.0 Intraoral X-ray units in November 2013 with a potential flaw in the mounting hardware.

Despite of our tight quality controls, we have found that a number of wall mount units were delivered with FOCUS™ X-ray machines that had not been properly welded at our supplier. As a result of the incomplete welding operation, the wall mount unit on the FOCUS™ X-ray unit may not properly support the unit as expected, which could cause the unit to drift. Depending on the extent of the shortened or missing welds, there is also a possibility that the wall mount unit will fail to support the weight of the unit, resulting in the FOCUS™ X-ray unit falling unexpectedly from its mounted position.

The serial numbers of the devices that potentially contain the wall mount units that were improperly welded fall between F21221 - F21643 (excluding F21642). The risk is limited to units in this serial number range.

As an authorized Instrumentarium dealer, we request that you carry out a field inspection on the suspected FOCUS™ 3.0 X-ray units to eliminate any potential risks related to the flawed wall mount units.

The suspected units with serial numbers between F21221 and F21643 (excluding F21642) need to be inspected for weld quality without any unnecessary delay. In cases where a defective wall mount unit is discovered, the FOCUS™ X-ray must be taken out of service and the wall mount unit replaced

INSTRUMENTARIUM

before continuing the use of the X-ray device. Replacement wall mount units will be available free of charge for the field repair. Please request service part 50864 when ordering the replacement part.

Please, refer to the attached technical bulletin (1-2014) for detailed instructions to inspect the FOCUS™ 3.0 X-ray unit wall mounts for the improper welds, as well as instructions to replace the wall mount unit in the event that a defective part is located.

As an authorized Instrumentarium Dental distributor we request you to take the following actions:

- Familiarize yourself with this document and its attachments.
- Return the Acknowledgment Letter confirming that you have received this letter and are planning to take appropriate actions to comply with the instructions above. Please return this document to Jussi Tanskanen at Instrumentarium Dental by the end of January, 2014
- Locate the suspected X-ray units in your region and organize an inspection of the units without a delay. Refer to Technical bulletin 1-2014 for details. If any defective units are discovered, replace the wall mount of the device under warranty before continuing the use of the FOCUS™ X-ray unit. Refer to this information letter in all communications to Instrumentarium Dental.
- Return the Confirmation Letter attached to this document, confirming that all of the FOCUS™ X-ray units between serial numbers F21221 and F21643 that have been delivered to your location have been inspected and/or corrected as required. Please return this form to Jussi Tanskanen at Instrumentarium Dental by the end of April, 2014.

We regret all the inconvenience caused by this issue. If there are any questions do not hesitate to contact us.

Mika Anttila

Marketing Director, General dentistry
Mika.anttila@instrudental.com
Tel. +358 40 768 3340

For more information you can contact also:

Tarja Päivinen

Director, Quality and Regulatory Affairs
Tarja.paivinen@instrudental.com
Tel. +358 40 719 0513

Mika Kyllönen

Aftersales Manager
Mika.Kyllonen@instrudental.com
Tel. + 358 40 080 5305

APPENDIXES

APPENDIX 1 Acknowledgement Letter
APPENDIX 2 Confirmation Letter

ACKNOWLEDGEMENT LETTER

We (distributor name) _____,

Confirm that we have read and understood the attached Field Safety Notice about FOCUS™ 3.0 field inspection and have informed the appropriate persons within our organization to ensure compliance with this request.

Date: _____

Distributor's representative signature: _____

Distributor's representative name : _____

Please, return a fulfilled acknowledgment letter back to Instrumentarium Dental by e-mail or mail.

RETURN ADDRESS:

Instrumentarium Dental
Jussi Tanskanen
Nahkelantie 160
FI-04300 FINLAND

Or

Jussi.tanskanen@instrudental.com

APPENDIX 2

CONFIRMATION LETTER

We, _____ (distributor name) confirm that

we have inspected the following FOCUS 3.0 intraoral X-ray units as required by the field safety notice attached to this confirmation letter.

We also confirm that actions have been taken as required by the field safety notice to replace any defective wall mount units that were identified through our inspection.

List of affected device serial numbers applicable:

Date: _____

Distributor's representative signature: _____

Distributor's representative name : _____

Please, return a fulfilled confirmation letter back to Instrumentarium Dental by e-mail or mail.

RETURN ADDRESS:

Instrumentarium Dental
Jussi Tanskanen
Nahkelantie 160
FI-04300 FINLAND

Or

Jussi.tanskanen@instrudental.com

| FSN n° 2014-14655 | |
|-------------------|-------------------------|
| File | 14655 |
| Manufacturer | Hospira Inc. |
| Device | GemStar infusion system |
| Commission date | 08/05/2014 |

15 January 2014

URGENT FIELD SAFETY NOTICE

Gemstar™ Infusion system (Versions prior to and including Phase 3)

| | |
|---------------|--|
| Product name: | GemStar 7 Therapy Pump – Versions prior to and including Phase 3 |
| | GemStar 6 Therapy Pump – Versions prior to and including Phase 3 |
| | GemStar Pain Management Pump – Versions prior to and including Phase 3 |
| List Number: | 13000, 13100, 13150 |
| EMEA FA ID: | Q.FA.EMEA.2013.038 |
| Date: | 15 th January 2014 |

Dear Healthcare Professional and Hospira Customer,

Hospira, Inc. (Hospira) is issuing this field safety notice to inform you about a situation that may occur during start up on GemStar (Versions prior to and including Phase 3) infusion pumps. **Please ensure that all potential users in your facility are made aware of this safety notification and the recommended actions.**

Issue: In GemStar (Versions prior to and including Phase 3) infusion pumps, list numbers 13000, 13100, and 13150, the connection between the beeper subassembly and the pump may fail. The GemStar infusion pump will identify this failure only during the “self-test” while powering up which will result in a Beeper Error (X100). All previous versions prior to Phase 3 carry the list numbers detailed within this FSN.

Risk to Health: If a GemStar fails the “self-test” upon power up, a delay of therapy may occur. The likelihood of harm will depend on the patient’s circumstances, such as the severity of their illness, the healthcare setting and the criticality of the infusion. When the GemStar is used in the home healthcare setting, there is no healthcare provider present, but it is extremely rare that a critical medication would be given in that setting. GemStar may also be used in the inpatient setting, where it is possible that a more critical infusion might be administered with this infusion pump; however, the healthcare professional to patient ratio is much higher for inpatients. The visual alarms on the pump make it likely that the delay in therapy will be noticed and of brief duration.

Affected Product Details: The products impacted by this issue are identified in the table below:

| Product Description | List Number |
|------------------------------|-------------|
| GemStar 7 Therapy Pump | 13000 |
| GemStar 6 Therapy Pump | 13100 |
| GemStar Pain Management Pump | 13150 |

Actions to be taken: There is no need to return your Gemstar devices at this time. Hospira recommends that users take the following actions:

1. Inform potential users in your organisation of this notification. Please also forward this notification to any organisation or persons (including home users) where the potentially affected devices have been transferred.
2. If a Beeper Error (X100) occurs during power up, remove the device from service immediately.
3. Complete the attached reply form and return it to the fax number or e-mail address on the form, even if you do not have the affected product.
4. If you have distributed the product further to the retail level, notify your accounts that received the product identified above of this notification and ask them to contact your local Hospira office to receive a reply form.

Product Correction: There is no further product correction at this time. As part of establishing a streamlined and modernised device portfolio that addresses customer needs, Hospira is in the process of retiring the GemStar family of infusion devices as communicated in the details of the global device strategy announced 1st May 2013. As of 31st July 2015, Hospira will consider the products within the GemStar Infusion System family retired and will no longer support them.

Please complete the attached Reply Form and return it to the fax number or e-mail address on the form, even if you do not have the affected product.

Hospira is committed to providing you with the highest level of service, product quality and reliability. We appreciate your understanding and we regret any inconvenience that may cause you.

Please forward this Field Safety Notice to all colleagues within your organisation who need to be aware of it or to any organisation or persons where the potentially affected devices have been transferred.

Please maintain awareness of this notice.

Should you have any further questions please do not hesitate to contact your local Hospira office:

| Hospira contact | Contact details | Areas of support |
|------------------------------------|--|---|
| Hospira EMEA Product Safety | T: +44 1926 834 400 Email to: devicecomplaintsemea@hospira.com | To report adverse events or product complaints |
| Hospira EMEA Quality | T: +31 36 5274 720 F: +31 36 5274 701 Email to: devicesfieldactions@hospira.com | Additional information and technical assistance |
| Local Contacts | | |

The Competent Authorities in all countries affected by this action have been informed of this field safety notice.

Yours sincerely,

Wilson Kennedy
EMEA Devices Quality Manager

URGENT FIELD NOTICE REPLY FORM

GemStar™ Infusion System (Versions prior to and including Phase 3)

| | |
|---------------|--|
| Product name: | GemStar 7 Therapy Pump GemStar 6 Therapy Pump GemStar Pain Management Pump |
| List Number: | 13000, 13100, 13150 |
| Hospira ref: | Q.FA.EMEA.2013.038 |

Section A

Hospital / Facility Details

Please fill out the information below and fax the completed form to Hospira at [local fax number].

| | |
|------------------------------|--|
| Name of Hospital / Facility: | |
| Hospital / Facility Address: | |
| Telephone Number: | |
| Name: | |
| Signature: | |
| Date: | |

Section B

☐ I have read and understood the contents of this Field Action, circulated it to all staff/departments that use this product and confirm that our inventory has been checked and we have no inventory of the listed products.

Or

Section C

☐ I have read and understood the contents of this Field Action, and circulated it to all staff/departments that use this product.

| FSN n° 2014-14726 | |
|-------------------|------------------|
| File | 14726 |
| Manufacturer | Synthes GmbH |
| Device | Softlock Forceps |
| Commission date | 08/05/2014 |

**To the ATTENTION of:
Operating Room Manager**

22 January 2014

URGENT MEDICAL DEVICE PRODUCT RECALL**Part Description / Part Number**

| Part Number | Part Description | Lot Numbers |
|--------------------|---|---------------------------------|
| 351.782 | Holding Forceps for SynReam Reaming Rod Ø 2.5 mm | 5923464, 5923181 |
| 399.091 | Bone Holding Forceps, self-centering, soft lock, length 191 mm | 5921007, 5921003, 5920633 |
| 399.092 | Reduction Forceps, large, with Points, soft lock, length 222 mm | 5923184, 5923086 |
| 399.097 | Bone Spreader, soft lock, width 8 mm, length 220 mm | 5923186 |
| 399.098 | Reduction Forceps, toothed, soft lock, length 194 mm | 5923278, 5923085 |

Dear Valued Customer,

Synthes GmbH is initiating a voluntary recall of the above mentioned articles and lots of the Softlock Forceps. Our records indicate that you may have inventory that is impacted by this recall.

Reason for the Recall:

The locking lever on the Softlock Forceps has the wrong angle and may not lock. There have been no reports of adverse events associated with this issue.

Potential hazard:

The presence of forceps that do not retain the desired position once applied may result in a moderate surgical delay, exposing the patient to the risk of prolonged exposure to anesthesia, the extension of an incision, the potential to delay healing, increase the risk for infection, and greater than anticipated postoperative pain. During use the surgeon may use forceps to retain an instrument, reduce a fracture, spread an osteotomy, or hold an implant in place. The soft lock feature present on forceps is utilized to retain a desired position at or near the target

operative area and if this feature fails the desired device, implant, or bone retention may be compromised. Causes of surgical delay as a result of this failure include the time taken to source an alternative instrument, and to reapply the new forceps. In the event that an alternative is not available, the surgeon may opt to utilize the malfunctioning forceps in a manual manner. This too may result in a surgical delay because the surgeon or other operative personnel may be required to apply manual tension to the device in order to maintain its desired position and achieve the operative goals. In addition, this may complicate the surgery because the surgeon may have to adjust his/her technique to accommodate the presence of additional hands near the operative field.

It is not reasonable to expect that the failure mode described will result in user injury or permanent impairment.

Customer immediate actions:

1. Immediately identify and quarantine all unused products listed above in a manner that ensures the affected products will not be used.
2. Review, complete, sign and return the attached reply form to your local DePuy Synthes sales organisation in accordance with the directions on the form within 5 business days of receipt of this notification.
3. Return any affected product within 30 business days. A credit note will be issued for the returned items.
4. Forward this notice to anyone in your facility that needs to be informed.
5. If any product listed below has been forwarded to another facility, contact that facility to arrange return.
6. Maintain awareness of this notice until all products listed below have been returned to DePuy Synthes.
7. Maintain a copy of this notice.

The applicable regulatory agencies are being notified. Synthes GmbH is voluntarily taking this action.

We apologise for any inconvenience that this product recall may create and appreciate your cooperation with our request. Should you have any queries please do not hesitate to contact your DePuy Synthes sales consultant.

Page 3 of 4

Thank you for your attention and cooperation.

Synthes GmbH

Claudia Allemann
Field Action Manager

Markus Wien
Director Quality Assurance Operations

Cc:

NOTICE: MEDICAL DEVICE RECALL R2013539**Softlock Forceps****Verification Section**

| Part Number | Part Description | Lot Numbers |
|--------------------|---|---------------------------------|
| 351.782 | Holding Forceps for SynReam Reaming Rod Ø 2.5 mm | 5923464, 5923181 |
| 399.091 | Bone Holding Forceps, self-centering, soft lock, length 191 mm | 5921007, 5921003, 5920633 |
| 399.092 | Reduction Forceps, large, with Points, soft lock, length 222 mm | 5923184, 5923086 |
| 399.097 | Bone Spreader, soft lock, width 8 mm, length 220 mm | 5923186 |
| 399.098 | Reduction Forceps, toothed, soft lock, length 194 mm | 5923278, 5923085 |

- ☐ We have located the identified product in stock; returned quantity is documented below, and a copy of this letter is being retained for our records.
- ☐ We do not have any identified product in stock; returned quantity is zero. We have retained a copy of this letter for our records.

RETURNED DEVICES (including quantity) and/or COMMENTS:

Hospital name: _____

Name/Title (please print) _____

Phone Number: _____

Signature and Date: _____

| | |
|-------------------|---|
| FSN n° 2014-14731 | |
| File | 14731 |
| Manufacturer | Micro-Tech Nanjing Co. Ltd. |
| Device | Single Use Biopsy Forceps/Einweg-Biopsiezange |
| Commission date | 08/05/2014 |

Urgent Field Safety Notice

Commercial name of the affected product: Single-use Biopsy Forceps

FSCA-identifier (e.g. date)

Type of action (e.g. chapter 4 definition of a FSCA).

Date: 2013.11.25

Attention: //////////////////////////////////

Details on affected devices:

The affected forceps products are come from Micro-Tech (Nanjing) Co., Ltd. The name is Single-use Biopsy Forceps and the type is NBF03-11123230(B130220001/B130107002/B130117001), NBF01-11123180(B130116002)

Description of the problem:

Include, as appropriate:

method of recovery, disposal or modification of device:

The endoscope was damaged and had to be repaired

The reason for FSN is that Micro-Tech received an accident on August 6th, 2013. The first complaint occurred in Germany on July 10th, 2013, and the description of it is: after taking the 4th biopsy probe the spike was sticking out (bended). The forceps was taken out with the endoscope and cut off. The endoscope was damaged and had to be repaired. The using product's REF at that time was NBF03-11023230, the LOT was B130220001. Our preliminary analysis is the material of needle is not stiffness. So products are all existed risk which the type is NBF03-11123230, NBF03-11123180, NBF01-11123230, NBF01-11123180, LOT is BXXXXXXXXX. Please check the products in store, if find the same type and LOT's products, please isolate and mark immediately, and mustn't sell. At the same time check the products that have been sold, whether have the existence of unused products. If find the same type and LOT's products, notify the customers to stop using immediately and help customers returns. Micro-Tech will arrange to recall the type and LOT products of, and Micro-Tech undertakes all costs that recall incurred.

Transmission of this Field Safety Notice: (if appropriate)

The notice needs to be passed to all distributors which sell this type and LOT in the market of EEA member states, Switzerland, Turkey and South Africa.

Contact reference person:

Name: **Adrian Zhang** Organisation: **Micro-Tech (Nanjing) CO., Ltd.**

Address: **NO.10 Gaoke Third Road, Nanjing National Hi-Tech, Industrial Development Zone, Nanjing 210061, Jiangsu Province, People's Republic of China**

Contact details: Tel: +86-25-58609879 Mail: adrian@micro-tech.com.cn

The undersign confirms that this notice has been notified the appropriate Regulatory Agency
(Closing paragraph)

Signature: Zhang

| | |
|-------------------|--------------------------|
| FSN n° 2014-14741 | |
| File | 14741 |
| Manufacturer | Siemens Medical Solution |
| Device | ACUSON SC2000 |
| Commission date | 08/05/2014 |

SIEMENS

CUSTOMER SAFETY ADVISORY NOTIFICATION

To users of the ACUSON SC2000™ ultrasound systems

Dear Valued Customer:

This letter is being sent to notify you of a potential malfunction with your ACUSON SC2000 ultrasound system with system serial numbers lower than 401100. If not addressed, this issue could lead to an unsafe mechanical condition when using your ACUSON SC2000 system.

When does this issue arise and what is the potential risk?

The user interface module on your ultrasound system consists of the console (monitor, control panel, and keyboard) and the swivel mechanism. The swivel mechanism is used for adjusting the height and rotation of the console and is attached to the column on your ultrasound system. If the connection of the user interface module to the column on the ultrasound system is loose, the entire module may fall off the ultrasound system.

In the event the module falls off the ultrasound system, there is potential risk of striking and harming anyone in the environment of the ultrasound system.

What steps can the user take to avoid potential risk of this issue?

If you observe looseness in the attachment of the user interface module to the column on the ultrasound system, immediately stop using your system and contact your Siemens service representative for repair.

How will this issue be resolved?


Your local Customer Service Engineer will contact you to make the necessary repairs on your system. You will be contacted to schedule an appointment. If you have any questions, please contact your local service support person for information regarding timelines and status.

Until the necessary repairs have been completed on your system, please share this information with all personnel within your organization who need to be aware of this issue.

As always, patient safety issues are a high priority, and we thus alert customers of product concerns which may occur under various conditions. To date, no patient injury has been reported. This problem was discovered as part of our ongoing quality process.

We sincerely regret any inconvenience this condition may cause in your daily operations.

Sincerely,

Sheila Pickering Ph.D.
Senior Director, Quality and Regulatory Affairs
Siemens Medical Solutions USA, Inc.
Ultrasound Business Unit

| | |
|-------------------|-------------------|
| FSN n° 2014-14851 | |
| File | 14851 |
| Manufacturer | Richard Wolf GmbH |
| Device | RIWOmobil |
| Commission date | 08/05/2014 |

Department: Technical Service

Clerk: Marek Rast

Tel.: +49 70 43 35 1046

Fax: +49 70 43 35 4467

Email: marek.rast@richard-wolf.com

Your ref. • Your communication dated • Our ref. VKT/ras • Date 17.02.2014

Urgent Safety Information

Safety measure relating to the RIWOmobil Product nos. 32116004, 32116005, 32116014 and 32116015

Addressed to:

Users, operators, distributors

Identification of the affected medical devices:

Equipment Cart **RIWOmobil product nos. 32116004, 32116005, 32116014 and 32116015.**
See enclosed list for affected series numbers.

Description of the problem:

Ladies and Gentlemen,

on the basis of a customer complaint, it has come to our notice that in two cases the power cable has been disconnected from the panel-mounted socket due to mechanical overload. This led to complete failure of the **RIWOmobil product no. 32116** equipment cart, as a result of the associated interruption of the power supply.

Measures to be taken by the addressees:

The product can continue to be used. However, we expressly draw your attention to the fact that until receipt and assembly of the new **Mounting Kit for Power Cables Product no. 643301191** the correct plug connection for the power cable to the panel-mounted socket should be checked before the equipment cart is used.

We hereby draw your attention to the fact that we have revised our Instruction Manual in light of this information. An updated version will be sent to you together with the mounting kit. You should replace the existing manual with the revised Instruction Manual.

In order to comply with the requirement to have a record of verification for the receipt of this **Urgent Safety Information**, we require the enclosed response form, and we ask you to return this to us by fax to the **fax number +49 7043 351360** by **04.03.2014**.

You are requested to retain this information at least until the measure has been carried out and completed.

Forwarding the information described here:

We ask you to ensure within your organization that all users of the products specified above and any other persons in a position to communicate relevant information are immediately informed about the content of this notice of Urgent Safety Information. If you have passed any products on to third parties, you are requested to forward a copy of this information notice or inform the contact persons listed below.

Measures by the manufacturer:

This Urgent Safety Information will be sent to all the affected customers by 18.02.2014.

Richard Wolf is modifying the equipment cart **RIWOmobil product no. 32116004, 32116005, 32116014 und 32116015**. We will send you the mounting kit and the revised Instruction Manual free of charge within the next 4 weeks including Mounting Instructions for the Power Cable Product no. 643301191.

Customer Service of instrument maker Richard Wolf will get in contact with you within 3 weeks of sending this letter and will then make arrangements for dispatch.

The responsible national government agency (Federal Institute for Drugs and Medical Devices in Bonn) has been informed about this **Urgent Safety Information**.

Your contacts:

For safety-related questions:

Mr. Heiko Sailer
Head of Quality Management and Regulatory Affairs
Tel.: +49 7043 35 1108
Fax: +49 7043 93 0048
Email: heiko.sailer@richard-wolf.com

For technical questions:

Mr. Marek Rast
Head of Technical Service
Tel.: +49 7043 35 1046
Fax: +49 7043 35 4467
Email: marek.rast@richard-wolf.com

We very much regret the inconvenience that you incur or have incurred as a result of the circumstances described above.

Yours sincerely,
Richard Wolf GmbH

p.p. Heiko Sailer
Head of Quality Management and Regulatory Affairs

p.p. Marek Rast
Head of Technical Service

Enclosure

| | |
|--------------------------|---|
| FSN n° 2014-14917 | |
| File | 14917 |
| Manufacturer | Trivascular Inc |
| Device | Ovation Prime Abdominal Stent Graft System |
| Commission date | 08/05/2014 |

Important Medical Device Information
Updated Product Instructions for Use (IFU)

for the Ovation Prime™ Abdominal Stent Graft System

February 27, 2014

Dear Clinical Partner,

As part of our commitment to continuous improvement, TriVascular is sending this communication to physician users of the Ovation Prime Abdominal Stent Graft System to provide notification of an update that will appear in a future revision of the product Instructions for Use (IFU). The IFU will be available in your country when the required translations and product registration/notification have been completed. This IFU update does not require rework or return of TriVascular product.

Details on the updated IFU:

TriVascular is incorporating an update to an existing IFU Caution statement in the Deployment Completion Instructions for the optional usage of angioplasty balloons to read as follows (edits indicated):

Caution: It is not recommended to balloon prior to 15 20 minutes after completion of the final polymer mix. Ballooning prior to 15 20 minutes could damage the sealing rings.

The change to the current ballooning time (from 15 to 20 minutes) is TriVascular's proactive response to two (2) reports, both in the United States of America, of damage to the sealing rings where ballooning occurred after 15 minutes and before 20 minutes. The two (2) patients experienced transient hypotension that stabilized after treatment in accordance with the current version of the product IFU, and the aneurysms for both patients were successfully excluded. These events represent a very low estimated world-wide incidence rate of less than 0.1%.

The change to ballooning time allows for additional curing of the fill polymer in variable clinical conditions. For example, if a patient has a low core body temperature (below 35°C) an additional minute of cure time per degree below 35°C may be required.

TriVascular commitment:

In our ongoing efforts to provide product education and guidance to physicians and also reduce potential patient safety risks, TriVascular updated the ballooning recommendations in the IFU.

TriVascular is dedicated to dynamically updating worldwide product instructions when new information is received from our Clinical Partners. We appreciate your willingness to work with us and are grateful to partner with you in the care of your patients. If you need additional information, please do not hesitate to contact your local TriVascular representative.

Sincerely,

Shari O'Quinn
VP of Clinical Affairs/Regulatory/Quality
TriVascular, Inc.

27 Feb 2014

URGENT FIELD SAFETY NOTICE - EXTENSION

| COMMERCIAL NAME OF AFFECTED PRODUCTS: | MATERIAL NUMBER | LOT |
|---------------------------------------|--|----------------------|
| Linton-Nachlas Tube | 204700 | Ref Appendix 2 |
| Sengstaken Tube | 204800 | |
| Pediatric Sengstaken Tube | 204801 | |
| Sengstaken-Blakemore Tube | 204802, 204803 | |
| Subcutaneous ureter bypass set | 335685 | |
| Endoguide | 680141, 680143, 680145, 680151, 680153, 680155 | |
| TYPE OF ACTION: | Recall (Extension) | |
| TELEFLEX REFERENCE: | 2013-08 | |

Dear Customer,

1. Details of affected devices

Teleflex has initiated a voluntary Field Safety Corrective Action for the above listed products.

2. Description of the problem

Teleflex Medical has issued a voluntary recall for the above listed products as the packaging may be compromised, and therefore the sterility of the product cannot be guaranteed. If non-sterile products are used, there is a possible risk of infection.

3. FIELD SAFETY CORRECTIVE ACTION INSTRUCTIONS:

ADVICE ON ACTION TO BE TAKEN BY MEDICAL STAFF

- We request that you check your inventory for product within the scope of this field action. Users should cease use and distribution of stock of affected product and quarantine immediately.
- If you do not have stock in scope of this field action as referred to in above table then mark the according checkbox on the Acknowledgement form (Appendix 1) and return the form to the fax number or e-Mail-address mentioned there.
- If you have stock from the affected product as referred to in above table, mark the according checkbox on the Acknowledgement form (Appendix 1). Contact customer service by calling the phone number mentioned in section 6 who will issue you with a return number. Write this return number into the respective field in the Acknowledgement form.
- Complete 'Appendix 1' for all products in your possession and under control. Return this form immediately to Customer Service.
- Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

INSTRUCTION FOR DISTRIBUTORS OF AFFECTED PRODUCT

1. If you are a distributor, provide this field safety notice to all of your customers who have received product in scope of this Field Action. Your customer is then required to complete the acknowledgement form and return this to you.
2. As a Distributor you are required to confirm to Teleflex that you have completed the field activity outlined above. Upon completion of your actions, please forward the completed Acknowledgement Form to Customer Service.
3. Please be aware that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
4. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/TK area, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

4. Teleflex

Teleflex informs all customers, employees of Teleflex and distributors on this Field Action.

5. Transmission of this Field Safety Notice

This notice should be passed on all persons who need to be aware within your organization or to any organization where the potentially affected devices have been transferred. Please consider end users, clinicians, risk managers, supply chain/distribution centres etc. in the circulation of this notice.

Maintain awareness of this notice until all required actions have been completed in your organisation

6. Contact reference person

Should you require any further information or support concerning this issue, please contact:

Customer Service:

Contact: Shane Kenny

FAX: +353 (0)1 4370773

Telephone: +353 (0)90 6460869

E-mail: orders.intl@teleflex.com

Please be advised that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities to which Teleflex distribute directly will be notified by Teleflex.

Teleflex is committed to providing high quality, safe and effective products. We sincerely apologize for any inconvenience this action may cause your operations. If you have any other questions, feel free to contact your local sales representative or Customer Service.

For and on behalf of Teleflex,

Padraig Hegarty

Padraig Hegarty

Senior Director of Quality International.

FIELD SAFETY CORRECTIVE ACTION ACKNOWLEDGEMENT FORM

PRODUCT FIELD ACTION BY TELEFLEX

IMMEDIATE ATTENTION REQUIRED

RETURN COMPLETED FORM IMMEDIATELY TO:

FAX: +353 (0)1 4370773

E-mail: orders.intl@teleflex.com

☐ We confirm receipt of this FSN and completed the required actions contained therein. We confirm that our inventory does **NOT** include products affected by this Field Action.

☐ We confirm receipt of this FSN and completed the required actions contained therein. We confirm our inventory does include products affected by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned.

Return Authorisation No _____

PLEASE PRINT PRODUCT QUANTITY NUMBERS CLEARLY.

| | | |
|--|---|-----------------|
| COMMERCIAL NAME OF AFFECTED PRODUCTS: | Linton-Nachlas Tube / Sengstaken Tube / Pediatric Sengstaken Tube / Pediatric Sengstaken Tube / Sengstaken-Blakemore Tube / Subcutaneous ureter bypass set/ Endoguide (EXTENSION) | |
| PRODUCT NUMBER | LOT NUMBER | QUANTITY |
| | | |
| | | |
| | | |
| | | |
| | | |
| <ul style="list-style-type: none"> Include a copy of the completed Acknowledgement Form in the returns package with the returned units Ensure the RAN number is clearly visible on the returns package. Please label returns as "Field Action Returns" | | |

Complete this Acknowledgement form and return immediately by using the fax number or e-mail address above.

| | |
|---|--------------------|
| INSTITUTION NAME (EG NAME OF HOSPITAL, HEALTH CARE ORGANISATION) | |
| | |
| INSTITUTION ADDRESS | Phone / Fax |
| | |
| FORM COMPLETED BY: | TITLE/ROLE |
| PRINT NAME: _____ | |
| SIGNATURE: _____ | |
| DATE | |

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| FSN n° 2014-15019 | |
|-------------------|------------------------------------|
| File | 15019 |
| Manufacturer | CareFusion Switzerland |
| Device | Alaris SmartSite Needle-Free Valve |
| Commission date | 08/05/2014 |

Name
Address

URGENT FIELD SAFETY NOTICE

Product Name: **Alaris® SmartSite® Needle-Free Valve**
Product Reference: **2000E7D**
Lot Numbers: **Various lots with Polycarbonate Material**
FSCA Identifier: **RA-2014-01-01**
Date: **February 2014**
Type of Action: **Removal & Destruction of Affected Products**

ATTENTION: Clinical Personnel, Risk Managers, Biomedical Personnel**Description of the Problem**

CareFusion has identified a potential risk with various lot numbers of Alaris® SmartSite® Needle-Free Valve - 2000E7D made with Polycarbonate. The affected SmartSite connector lots may unintentionally disconnect from a female Luer, may be difficult to disconnect from a female Luer, or may fail to disconnect from a female Luer once attached. In rare occasions limited leakage may also be observed if the connector disconnects from the female Luer during infusion. An inability to disconnect may require replacement of a line.

Note: This recall does NOT involve MaxPlus or MaxZero connectors.

Potential Risk

A disconnection from the female Luer can cause a delay of infusion, an interruption of infusion, an under infusion or potentially air entering the fluid path and/or limited leakage. An inability to disconnect can cause delay in infusion or an interruption of infusion. We have not received any reports of serious injuries or death related to this issue.

Root Cause

The root cause is still under investigation. Preliminary analysis has indicated that SmartSite with polycarbonate material may cause incompatibility issues with the male Luer or thread housing when attaching to a female Luer.

SmartSite with acrylic material does not appear to cause these incompatibility issues and is not affected by this FSCA. There is a visible colour difference between the polycarbonate and acrylic base. The polycarbonate material is bluish-grey whilst the acrylic material is clear. As shown in the pictures below.

Polycarbonate Base

Acrylic Base

Products Potentially Affected

Lot numbers of SmartSite with polycarbonate are as follows:

13086202, 13086263, 13086309, 13095185, 13095913, 13095914, 13096487,
13096488, 13096489, 13105366, 13105478, 13106563, 13106710, 13106711,
13115316

Our traceability analysis determined that you have received SmartSite with polycarbonate material and therefore might experience this defect.

Action Required

In order that the potentially affected products are removed from use and destroyed please follow instructions below:

| Step | Action | |
|------|---|--|
| 1 | Inspect inventory for product affected lot numbers of 2000E7D. (listed in Appendix 1) | |
| | If... | Then... |
| | No affected lot numbers are found | <ul style="list-style-type: none"> Complete sections A & B of Appendix 1 |
| | Affected lot numbers are found | <ul style="list-style-type: none"> Complete sections A & C of Appendix 1 Destroy affected product according to hospital protocol |
| 2 | Return completed verification form (Appendix 1) to your CareFusion representative no later than 31 May 2014 | |
| 3 | Please contact your local CareFusion representative to order alternative products. | |

On receipt of the completed verification form Appendix 1, CareFusion will credit you for any product destroyed.

Your competent authority has already been notified of this Field Safety Corrective Action by CareFusion.

We sincerely apologise for any inconvenience this action may have caused you or your staff.

Should you have any questions or require assistance relating to this Field Safety Corrective Action, please contact your Local CareFusion representative.

Transmission of this Field Safety Notice

Please distribute this notice to all those who need to be aware of this action within your organisation.

Sincerely,

CareFusion Representative

Appendix 1

URGENT FIELD SAFETY NOTICE – Verification Form

Product Name: **Alaris® SmartSite® Needle-Free Valve**
 Product Reference: **2000E7D – Various Lots**
 FSCA Identifier: **RA-2014-01-01**
 Date: **February 2014**
 Type of Action: **Removal & Destruction of Affected Products**

Section A

| | |
|--|--|
| Name of Hospital / Facility | |
| Hospital / Facility Address | |
| Name of distributor (if applicable) | |
| Telephone Number | |
| Name | |
| Signature | |
| Date | |

Section B

☐ I have read and understood the contents of this Field Action and confirm that our inventory has been checked and we have no inventory of the listed products.

Section C

☐ I have read and understood the contents of this Field Action and confirm that our inventory has been checked and the following products have been destroyed:

| Batch Number | Quantity Destroyed | Batch Number | Quantity Destroyed | Batch Number | Quantity Destroyed |
|--------------|--------------------|--------------|--------------------|--------------|--------------------|
| 13086202 | | 13095914 | | 13105478 | |
| 13086263 | | 13096487 | | 13106563 | |
| 13086309 | | 13096488 | | 13106710 | |
| 13095185 | | 13096489 | | 13106711 | |
| 13095913 | | 13105366 | | 13115316 | |

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Nom du Centre :

Signature, date :

| | |
|--------------------------|--|
| FSN n° 2014-15103 | |
| File | 15103 |
| Manufacturer | Lutonix Inc. C.R. Bard |
| Device | Lutonix Drug Coated PTA Dilatation Catheter |
| Commission date | 08/05/2014 |

«Hospital_Name»

«Users_Name» - «Department»

«Customer_Address»

«Zip_Code» «City» - «Country_name»

Reference: 90926750-FA

March 24, 2014

Field Safety Notice (for Distributor) Urgent Medical Device Recall

Teleflex Medical CAPIO™ Sutures

Dear «Users_Name»,

Boston Scientific, as the distributor of **CAPIO™ Sutures**, is implementing this Medical Device Recall at the request of Teleflex Medical, the manufacturer and owner of these devices. Boston Scientific has been notified by Teleflex Medical that certain lots of their products did not meet minimum needle attachment strength, knot tensile strength, minimum diameter, or minimum resorption strength requirements. This may result in suture breakage, delay in procedure, wound dehiscence, or the need for additional surgical intervention.

Our records indicate that your facility received some of the concerned product. **The table below provides a complete list of all affected products**, including Product Description, Material Number (UPN), Lot/Batch numbers and Expiration date. **Please note that only the material listed in the table below is affected. No other Boston Scientific product is involved by this Field Safety Notice.**

Further distribution or use of any remaining product affected by this action should cease immediately.

| Product Description | Material Number (UPN) | | Catalog # | Lot | Expiration Date |
|---|-------------------------------|----------------------|-----------|--|------------------------------------|
| | Outer carton label of 12 pack | Inner Pouch (single) | | | |
| CAPIO™ Braided Polyester Green PTFE Impregnated Polyester Fiber, Nonabsorbable Surgical Suture | M0068331141 | M0068331140 | 833-114 | 02A1202112, 02C0902841, 02C0903374, 02D1202794 02F0902697, 02K1100404 02L1202369, 02M0802509 | January 1, 2012 to August 31, 2018 |
| CAPIO™ Polypropylene Blue Monofilament Nonabsorbable Surgical Sutures | M0068331231 | M0068331230 | 833-123 | 02A0902858, 02D0900775 02G1003092, 02H1100535 02H1100536, 02K1102309 | January 1, 2012 to August 31, 2018 |

| Product Description | Material Number (UPN) | | Catalog # | Lot | Expiration Date |
|---|-------------------------------|----------------------|-----------|--|--|
| | Outer carton label of 12 pack | Inner Pouch (single) | | | |
| CAPIO™ Polypropylene Blue Monofilament Nonabsorbable Surgical Sutures | M0068331241 | M0068331240 | 833-124 | 02A1200503, 02B1101450 02C0903400, 02D0900010 02D0900103, 02D1003285 02D1003286, 02F1302321 02G1000045, 02G1300348 02G1301122, 02H1100494 02L0803407, 02L1100009 02M1002289, 02M1101933 | January 1, 2012 to August 31, 2018 |
| CAPIO™ Bondek® Green Braided-Coated Polyglycolic Acid Sutures, Synthetic Absorbable Surgical Sutures | M0068331361 | M0068331360 | 833-136 | 02H1300693, 02L1200939 02L1200940 | January 1, 2012 to August 31, 2018 |
| CAPIO™ Monodek® Violet Monofilament Polydioxanone Sutures, Absorbable Surgical Sutures | M0068331371 | M0068331370 | 833-137 | 02A0901938, 02B0901537 02C1102979, 02C1102983 02C1302078, 02D1300294 02D1301164, 02D1302468 02H1001025, 02H1001026 02H1001027, 02H1103434 02H1200349, 02K1201354 | January 1, 2012 to August 31, 2018 |
| CAPIO™ Bondek® Green Braided-Coated Polyglycolic Acid Sutures, Synthetic Absorbable Surgical Sutures | M0068332131 | M0068332130 | 833-213 | 02H1100687 | January 1, 2012 to August 31, 2018 |

INSTRUCTIONS:

- 1. Check your Inventory for the products affected by this Field safety Notice.**
- 2. Please segregate and do not distribute any of the affected products found in your inventory**
- 3. Please notify all your customers that have received affected product of this Field Safety Notice.** To effectively manage this Field Safety Notice, your accounts are to communicate directly with you, not Boston Scientific. If any of your customers are distributors, please notify them that they must communicate this recall to the medical facility level.
- 4. Have all affected product returned to your facility and held in quarantine for reconciliation**
- 5. Please complete the attached Verification Form even if you do not have any product to return.**
- 6. When completed, please fax the Verification Form** to your local Boston Scientific Office to the attention of «Customer_Service_Fax_Number» **on or before April 7, 2014.**
- 7. If you have products to return, please package them in appropriate shipping box and contact «Customer_Service_Tel» of your local Boston Scientific Office, to arrange return.**

You will receive Replacement for all recalled product that is returned to Boston Scientific.

Your Competent Authority is being notified (by Teleflex Medical) of this Field Safety Notice.

We regret any inconvenience that this action may cause and we appreciate your understanding as we take action to ensure patient safety and customer satisfaction.

If you have any questions or would like assistance with this Field Safety Notice, please contact your local Sales Representative.

Yours sincerely,

Marie Pierre Barlanga
Quality Department
Boston Scientific International S.A.

Attachment: Verification Form

| | |
|--------------------------|--|
| FSN n° 2014-15160 | |
| File | 15160 |
| Manufacturer | Temena SAS |
| Device | Ultrasound needles, USB and USC |
| Commission date | 08/05/2014 |

Carrières sur Seine, 2014-03-24

**Urgent recall of Temena USB and USC needles,
article number: USB050-22 and USC050-22**

Dear Ladies and Gentlemen,

We got feedback from the market that some of our USB and USC needles were leaking at the hub. We cannot exclude that the some devices of the batches you received might be affected. Temena SAS therefore decided to recall the concerned medical devices voluntarily.


Following products/batches are affected:

| Product designation | Manufacturer | Article number | Batch number | Expiry date |
|--|--|----------------|--------------|-------------|
| Ultra sound needle 22G (0.7 mm) x 50 mm | Temena SAS, Carrières sur Seine, France | USC050-22 | B3170005 | 2018-05-31 |
| Ultra sound needle, 30° bevel | | USB050-22 | B3070012 | 2018-02-28 |
| Ultra sound needle, 30° bevel | | USB050-22 | B3250001 | 2018-06-30 |
| Ultra sound needle, 30° bevel | | USB050-22 | B3250026 | 2018-06-30 |

IMPORTANT INSTRUCTIONS:

- ☐ Ensure that the devices questions will not be applied to humans.
- ☐ Opened packages shall not be used any more.
- ☐ Please send us back the request of products (including opened packages)
- ☐ Inform your distribution partners / further resellers about this recall to achieve a 100% recall of all unused products.
- ☐ Please add a copy of our invoice / respectively of our delivery notes, to carry out a cost reimbursement or receive compensation delivery or a credit note.
- ☐ Please use attached form to answer Temena SAS.

We sincerely apologize for the inconvenience this may have caused by this product recall we thank you in advance for your support.

.....
Wilhelm Waskönig
President
.....
Dr. Bernhard Kusian
vigilance manager

Response Product Recall

Please reply to Fax Nr. +33 (0)1.30.86.05.31

via email: info@temena.fr

**Re: Urgent recall of Temena USB (REF USB050-22) and
USC needles (REF USC050-22)**

Affected products:

| Product designation | Manufacturer | Article number | Batch number | Expiry date |
|--|--|----------------|--------------|-------------|
| Ultra sound needle 22G (0.7 mm) x 50 mm | Temena SAS, Carrières sur Seine, France | USC050-22 | B3170005 | 2018-05-31 |
| Ultra sound needle, 30° bevel | | USB050-22 | B3070012 | 2018-02-28 |
| Ultra sound needle, 30° bevel | | USB050-22 | B3250001 | 2018-06-30 |
| Ultra sound needle, 30° bevel | | USB050-22 | B3250026 | 2018-06-30 |

Customer:

Company:
Street:
ZIP / Town:

| | | | | |
|--------------------------|--|------|---------|--|
| <input type="checkbox"/> | We have used all devices of the concerned LOTs. No defects have been observed so far. | | | |
| <input type="checkbox"/> | We still have following qualities of the concerned LOTs in our stock: | | | |
| <input type="checkbox"/> | REF: | LOT: | Amount: | |
| <input type="checkbox"/> | REF: | LOT: | Amount: | |
| <input type="checkbox"/> | REF: | LOT: | Amount: | |
| <input type="checkbox"/> | REF: | LOT: | Amount: | |
| <input type="checkbox"/> | We will return our remaining products from these LOTs, including opened cardboxes. | | | |
| <input type="checkbox"/> | | | | |
| <input type="checkbox"/> | | | | |

Please check where appropriate.

Please send us

a credit note ☐, a cost reimbursement ☐, a product exchange ☐.

Further comments:

You will reach me as signatory under the following phone number.

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SUNRISE MEDICAL LOGISTICS B.V.